GenoType[®] MTBDRsl

VER 1.0

CE

Deutsch: S. 2-7

English: p. 8-14

GenoType® MTBDRsI

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GenoType® MTBDRsl

Molecular Genetic Assay for Identification of Resistances to Fluoroquinolones, Aminoglycosides/Cyclic Peptides and Ethambutol of the *Mycobacterium tuberculosis* Complex

Methodology

The GenoType® MTBDRsl test is based on the DNA•STRIP® technology and permits the molecular genetic identification of the *M. tuberculosis* complex and its resistance to fluoroquinolones (e.g. ofloxacin and moxifloxacin) and/or aminoglycosides/cyclic peptides (injectable antibiotics as capreomycin, viomycin/kanamycin, amikacin) and/or ethambutol from cultivated samples or pulmonary smear-positive clinical specimens. The identification of resistance to fluoroquinolones is enabled by the detection of the most significant mutations of the *gyrA* gene (coding for DNA gyrase). For detection of resistance to aminoglycosides/cyclic peptides, the 16S rRNA gene (*rrs*) and for detection of resistance to ethambutol the *embB* gene (which, together with the genes *embA* and *embC*, codes for arabinosyl transferase) are examined.

The whole procedure is divided into three steps: DNA extraction from cultured material (solid/liquid medium) or clinical specimens (pulmonary, smear-positive, decontaminated) – the necessary reagents are not provided, a multiplex amplification with biotinylated primers (the necessary thermostable DNA polymerase is not provided), and a reverse hybridization.

The hybridization includes the following steps: chemical denaturation of the amplification products, hybridization of the single-stranded, biotin-labeled amplicons to membrane-bound probes, stringent washing, addition of a streptavidin/alkaline phosphatase (AP) conjugate, and an AP mediated staining reaction. A template ensures the easy and fast interpretation of the banding pattern obtained.

Storage and Precautions

Store Primer/Nucleotide Mix (PNM) at 2-8°C upon arrival isolated from any potential source of contaminating DNA. If longer storage (more than 4 weeks) is required, store at -20°C. In order to avoid repeated freezing and thawing, aliquot PNM. Store all other kit components at 2-8°C. Do not use the reagents beyond their expiry date.

Patient specimens and cultures made from patient specimens must always be considered as potentially infectious. Samples from risk patients and cultures made from those samples must always be labeled and handled under suitable safety conditions. Observe all federal, state, and local safety and environmental regulations. Always wear suitable protective clothing and gloves. Specimen treatment and sample preparation up to and including the heat inactivation step must be carried out in a class II safety cabinet. Before the heat inactivation step, samples must be centrifuged in an aerosol-tight rotor. Open aerosol-tight rotor in safety cabinet only. After heat inactivation, a standard rotor can be used for spinning the samples outside the safety cabinet. Observe the usual precautions for amplification set-up. It is essential that all reagents and materials used for DNA extraction and amplification set-up are free from DNases.

When handling kit reagents, the following special safety measures must be applied:

The **Denaturation Solution** (DEN) contains <2% NaOH and is irritating to eyes and skin (R 36/38 and S 26-37/39-45).

The **Substrate Concentrate** (SUB-C) contains Dimethyl Sulfoxide and is irritating (R 36/37/38, S 23-26-36).

For additional information, please refer to material safety data sheets which can be downloaded from: www.hain-lifescience.com/products/msds.html

Quality Control

In order to validate the correct performance of the test and the proper functioning of reagents, each strip includes 5 control zones:

- a Conjugate Control zone to check the binding of the conjugate on the strip and a correct chromogenic reaction
- an Amplification Control zone to check for a successful amplification reaction
- three Locus Control zones (*gyrA*, *rrs* and *embB*) checking the optimal sensitivity of the reaction for each of the tested gene loci

DNA Extraction

Bacteria grown on solid medium (e. g. Loewenstein-Jensen, Middlebrook) or in liquid medium (e. g. BACTEC, MB-Check) may be used, as well as smear-positive clinical specimens (pulmonary samples). The test must not be used to detect mycobacteria directly from smear-negative patient specimens. The working area must be free from amplified DNA. It is crucial to heat samples to 95°C for at least 20 min in order to inactivate vegetative bacteria. Any DNA extraction procedure producing amplifiable DNA from bacteria can be used. The following quick protocol normally yields DNA suitable for amplification:

- 1a. When using bacteria grown on solid medium, collect bacteria with an inoculation loop and suspend in approximately 300 µl of water (molecular biology grade).
- 1b. When using bacteria grown in liquid media, directly apply 1 ml, when using clinical specimens, apply 500 µl of a decontaminated* sample. Pellet bacteria by spinning for 15 min in a standard table top centrifuge with an aerosol-tight rotor in a class II safety cabinet at approximately 10000 x g. Discard supernatant and resuspend bacteria in 100-300 µl of water (for culture samples), or 100 µl of water (for clinical specimens) by vortexing.
2. Incubate bacteria from 1a or 1b for 20 min at 95°C in a water bath.
3. Incubate for 15 min in an ultrasonic bath (for cultured material this step is optional).
4. Spin down for 5 min at full speed and directly use 5 µl of the supernatant for PCR. In case DNA solution is to be stored for an extended time period, transfer supernatant to a new tube.

* Samples must be processed using the NALC/NaOH method according to the CDC publication "Public Health Mycobacteriology: A Guide for the Level III Laboratory".

Detailed protocols can be requested from your local distributor or through: www.hain-lifescience.com

Amplification

Prepare the amplification mix (45 µl) in a DNA-free room. The DNA sample should be added in a separate area.

Per tube mix:

- 35 µl PNM – provided
- 5 µl 10x polymerase incubation buffer – not provided
- x µl MgCl₂ solution¹¹ – not provided
- 1-2 unit(s) thermostable DNA polymerase (refer to manual) – not provided
- y µl water to obtain a volume of 45 µl (not considering volume of enzyme) – not provided
- Add 5 µl DNA solution leading to a final volume of 50 µl (not considering volume of enzyme).

¹¹ Depending on the enzyme/buffer system used, the optimal MgCl₂ concentration may vary between 1.5 and 2.5 mM. Please note that some incubation buffers already contain MgCl₂.

The performance evaluation of the **GenoType® MTBDRsI** assay was carried out using the HotStarTaq DNA Polymerase from Qiagen. When using this enzyme, the following amounts are required per sample:

- 35 µl PNM – provided
- 5 µl 10x PCR Buffer for HotStarTaq (contains 15 mM MgCl₂) – not provided
- 2 µl 25 mM MgCl₂ solution – not provided
- 0.2 µl (1 U) HotStarTaq – not provided
- 3 µl water (molecular biology grade) – not provided
- 5 µl DNA solution (add in a separated area)

The final MgCl₂ concentration in this amplification mix is 2.5 mM.

Determine the number of samples to be amplified (number of samples to be analyzed plus control samples). A negative control sample, for example, contains 5 µl of water instead of DNA solution. Prepare a master mix containing all reagents except for DNA solution and mix well (do not vortex). Aliquot 45 µl in each of the prepared PCR tubes.

Amplification profile²¹:

		culture samples	clinical specimens
15 min	95°C	1 cycle	1 cycle
30 sec	95°C	} 10 cycles	10 cycles
2 min	58°C		
25 sec	95°C	} 20 cycles	30 cycles
40 sec	53°C		
40 sec	70°C		
8 min	70°C	1 cycle	1 cycle

²¹ Applies to the Taq polymerase used for validation. When using other hot start DNA polymerases, the time interval of the first step might have to be reduced (please refer to manual of the enzyme).

Amplification products can be stored at +4 to -20°C.

For checking the amplification reaction, 5 µl of each sample might be directly applied to a 2% agarose gel without the addition of loading buffer. The amplicons have a length of approximately 63 bp (Amplification Control), 115 bp (*M. tuberculosis* complex), 123 bp (*gyrA*), 263 bp (*rrs*), and 138 bp (*embB*), respectively.

Hybridization

Preparation

Prewarm shaking water bath/TwinCubator® to 45°C; the maximum tolerated deviation from the target temperature is +/-1°C. Prewarm solutions HYB and STR to 37-45°C before use. The reagents must be free from precipitates (note, however, that solution CON-D is opaque). Mix if necessary. Warm the remaining reagents with the exception of CON-C and SUB-C to room temperature. Using a suitable tube, dilute Conjugate Concentrate (CON-C, orange) and Substrate Concentrate (SUB-C, yellow) 1:100 with the respective buffer (**CON-C with CON-D, SUB-C with SUB-D**) in the amounts needed. Mix well and bring to room temperature. For each strip, add 10 µl concentrate to 1 ml of the respective buffer. Dilute CON-C before each use. Diluted SUB-C is stable for 4 weeks if stored at room temperature and protected from light.

1. **Dispense 20 µl of Denaturation Solution (DEN, blue) in a corner of each of the wells used.**

2. **Add to the solution 20 µl of amplified sample, pipette up and down to mix well and incubate at room temperature for 5 minutes.**

Meanwhile, take strips out of the tube using tweezers and mark them with a pencil underneath the colored marker. Always wear gloves when handling strips.

3. **Carefully add to each well 1 ml of prewarmed Hybridization Buffer (HYB, green). Gently shake the tray until the solution has a homogenous color.**

Take care not to spill solution into the neighboring wells.

4. **Place a strip in each well.**

The strips must be completely covered by the solution and the coated side (identifiable by the colored marker near the lower end) must face upward. Using tweezers, turn over strips which might have turned when immersed in the solution. Carefully clean tweezers after each use to avoid contamination. This also applies to all following steps.

5. **Place tray in shaking water bath/TwinCubator® and incubate for 30 minutes at 45°C.**

Adjust the shaking frequency of the water bath to achieve a constant and thorough mixing of the solution. To allow adequate heat transfer, the tray must be dipped into the water to at least 1/3 of its height.

6. **Completely aspirate Hybridization Buffer.**

For example, use a Pasteur pipette connected to a vacuum pump.

7. **Add 1 ml of Stringent Wash Solution (STR, red) to each strip and incubate for 15 minutes at 45°C in shaking water bath/TwinCubator®.**

8. **Work at room temperature from this step forward.**

Completely remove Stringent Wash Solution.

Pour out Wash Solution in a waste container and remove all remaining fluid by turning the tray upside down and gently striking it on an absorbent paper.

This also applies to all other wash steps.

9. **Wash each strip once with 1 ml of Rinse Solution (RIN) for 1 minute on shaking platform/TwinCubator® (pour out RIN after incubation).**

10. **Add 1 ml of diluted Conjugate (see above) to each strip and incubate for 30 minutes on shaking platform/TwinCubator®.**

11. **Remove solution and wash each strip twice for 1 minute with 1 ml of Rinse Solution (RIN) and once for 1 minute with approx. 1 ml of distilled water (e. g. use wash bottle) on shaking platform/TwinCubator® (pour out solution each time).**

Make sure to remove any trace of water after the last wash.

12. **Add 1 ml of diluted substrate (see above) to each strip and incubate protected from light without shaking.**


Depending on the test conditions (e. g. room temperature), the substrate incubation time can vary between 3 and 20 minutes. Extended substrate incubation times can lead to increased background staining and might impair interpretation of the results.

13. **Stop reaction by briefly rinsing twice with distilled water.**

14. **Using tweezers, remove strips from the tray and dry them between two layers of absorbent paper.**

Evaluation and Interpretation of Results

Paste strips and store protected from light. An evaluation sheet is provided with the kit. When using this evaluation sheet, paste the developed strips in the designated fields by aligning the bands CC and AC with the respective lines on the sheet. Determine the resistance status and note down in the respective column; as an help for interpretation, evaluation examples are given in the subsequent chapter. The supplied template also serves as an aid for evaluation and must be aligned with the bands CC and AC of the strip as well. Each strip has a total of 22 reaction zones (see figure).

Note: The strip is not displayed in original size.

Not all bands of a strip have to show the same signal strength.

Conjugate Control (CC)

A line must develop in this zone, documenting the efficiency of conjugate binding and substrate reaction.

Amplification Control (AC)

When the test is performed correctly, a control amplicon will bind to the Amplification Control zone. If this band is developed, mistakes during extraction and amplification setup and the carry-over of amplification inhibitors can be excluded.

In case of a positive test result, the signal of the Amplification Control zone can be weak or even vanish totally. This might be due to competition of the single reactions during amplification. In this case, however, the amplification reaction was performed correctly and the test does not have to be repeated.

A missing AC band in case of a negative test result indicates mistakes during amplification set-up, or carry-over of amplification inhibitors. In this case, the test is not valid and the respective sample has to be repeated.

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This zone hybridizes, as known, with amplicons generated from all members of the *Mycobacterium tuberculosis* complex. If the TUB zone is negative, the tested bacterium does not belong to the *M. tuberculosis* complex and cannot be evaluated by this test system.

Locus Controls (*gyrA*, *rrs* and *embB*)

The Locus Control zones detect a gene region specific for the respective locus and must always stain positive. If neither the Locus Control probe nor the Wild type or Mutation probes of one of the three genes examined are developed, the test can not be evaluated.

Wild type probes

The wild type probes comprise the most important resistance regions of the respective genes [see figure 1 as well as tables 1, 2, and 3]. When all wild type probes of a gene stain positive, there is no detectable mutation within the examined regions. Hence, the strain tested is sensitive for the respective antibiotics.

In case of a mutation, the respective amplicon cannot bind to the corresponding wild type probe. The absence of a signal for at least one of the wild type probes hence indicates a resistance of the tested strain to the respective antibiotics.

Only those bands whose intensities are about as strong as or stronger than that of the Amplification Control zone are to be considered.

Each pattern deviating from the wild type pattern indicates a resistance of the tested strain. The banding pattern obtained with the *gyrA* probes allows to draw a conclusion about a resistance to fluoroquinolones (e.g. ofloxacin or moxifloxacin), the banding pattern obtained with the *rrs* probes allows to draw a conclusion about a resistance to aminoglycosides (e.g. capreomycin or viomycin)/cyclic peptides (e.g. kanamycin or amikacin) and the banding pattern obtained with the *embB* probes allows to draw a conclusion about a resistance to ethambutol of the strain tested.


Fig. 1: Fluoroquinolone resistance region of the *gyrA* gene

gyrA WT1-3: *gyrA* wild type probes; *gyrA* MUT1-3: *gyrA* mutation probes. The numbers specify the positions of the amino acids (codons) for all mutations listed in the table. The codons for which mutation probes were designed are highlighted.

Mutation probes

The mutation probes detect some of the most common resistance mediating mutations [see table 1, 2, and 3].

Only those bands whose intensities are about as strong as or stronger than that of the Amplification Control zone are to be considered.

Each pattern deviating from the wild type pattern indicates resistance of the tested strain. The banding pattern obtained with the *gyrA* probes allows to draw a conclusion about a resistance to fluoroquinolones (e.g. ofloxacin or moxifloxacin), the banding pattern obtained with the *rrs* probes allows to draw a conclusion about a resistance to aminoglycosides (e.g. capreomycin or viomycin)/cyclic peptides (e.g. kanamycin or amikacin) and the banding pattern obtained with the *embB* probes allows to draw a conclusion about a resistance to ethambutol of the strain tested.

Note the following special cases:

1. There is a possibility that the specimen tested contains a heteroresistant strain. In case of a heteroresistance, a mutated as well as a wild-type sequence can be detected in the respective sample; hence one of the mutation probes as well as the corresponding wild type probe may stain positive on the respective strip. Whether the respective resistance becomes phenotypically evident depends on the ratio of mutated and non-mutated sequences at investigation.
2. There is a possibility that the tested specimen contains more than one *M. tuberculosis* complex strain [due to mixed culture or contamination]. If at least one of these strains harbors a mutation, one of the mutation probes as well as the corresponding wild type probe may stain positive. Whether the respective resistance becomes phenotypically evident, depends on the ratio of resistant and sensitive strain at investigation.

Table 1: Mutations in the gene *gyrA* and the corresponding wild type and mutation probes

(according to Cheng *et al.* 2004, Aubry *et al.* 2006, Matrat *et al.* 2006 and Kocagöz *et al.* 1996)

failing wild type probe(s)	codons analyzed	mutation probe	mutation	phenotypic resistance* ¹
<i>gyrA</i> WT1	85-90		C88S	FLQ ^R
<i>gyrA</i> WT1	85-90		A88T	FLQ ^R
<i>gyrA</i> WT2	89-93	<i>gyrA</i> MUT1	A90V	FLQ ^R
<i>gyrA</i> WT2	89-93	<i>gyrA</i> MUT2	S91P	FLQ ^R
<i>gyrA</i> WT3	92-97	<i>gyrA</i> MUT3A	D94A	FLQ ^R
<i>gyrA</i> WT3	92-97	<i>gyrA</i> MUT3B	D94N	FLQ ^R
<i>gyrA</i> WT3	92-97	<i>gyrA</i> MUT3B	D94Y	FLQ ^R
<i>gyrA</i> WT3	92-97	<i>gyrA</i> MUT3C	D94G	FLQ ^R
<i>gyrA</i> WT3	92-97	<i>gyrA</i> MUT3D	D94H* ²	FLQ ^R

*¹ FLQ: fluoroquinolones (e.g. ofloxacin and moxifloxacin); ^R = resistance

*² This rare mutation has only been detected theoretically (*in silico*) yet. It is therefore possible that it cannot be detected *in vitro*.

Table 2: Mutations in the gene *rrs* and the corresponding wild type and mutation probes

(according to Maus *et al.* (Feb) 2005 and Maus *et al.* (Aug) 2005)

failing wild type probe(s)	codons analyzed	mutation probe	mutation	phenotypic resistance* ³
<i>rrs</i> WT1	1401	<i>rrs</i> MUT1	A1401G	CAP ^R , Vio ^S , AMK ^R , KAN ^R
<i>rrs</i> WT1	1402	-	C1402T	CAP ^R , Vio ^R , AMK ^S , KAN ^R
<i>rrs</i> WT2	1484	<i>rrs</i> MUT2	G1484T	CAP ^R , Vio ^R , KAN ^R , AMK ^R

*³ CAP: capreomycin, Vio: viomycin, KAN: kanamycin, AMK: amikacin

^R = resistance; ^S = sensitivity

Table 3: Mutations in the gene *embB* and the corresponding wild type and mutation probes (according to Johnson *et al.* 2006 and Plinke *et al.* 2006)

failing wild type probe(s)	codons analyzed	mutation probe	mutation	phenotypic resistance ^{*4}
<i>embB</i> WT	306	<i>embB</i> MUT1A	M306I ^{*5}	EMB ^R
<i>embB</i> WT	306	<i>embB</i> MUT1B	M306V	EMB ^R
<i>embB</i> WT	306	-	M306I ^{*6}	EMB ^R
<i>embB</i> WT	306	-	M306I ^{*7}	EMB ^R

^{*4} EMB: ethambutol

^{*5} M306I: base exchange at codon 306: ATG→ATA

^{*6} M306I: base exchange at codon 306: ATG→ATC

^{*7} M306I: base exchange at codon 306: ATG→ATT

^R = resistance

Evaluation Examples


Figure 2: Examples for banding patterns and their evaluation with respect to resistances to fluoroquinolones and/or aminoglycosides/cyclic peptides and/or ethambutol

If all wild type bands display a signal, this is classified as positive and marked in the WT column of the respective gene as "+". If at least one of the wild type bands is absent, this is classified as negative and marked in the WT column as "-". Negative entries are only made to the mutation columns when none of the mutation bands display a coloration. If at least one of the mutation bands display a coloration, this is classified as positive and the MUT column of the respective gene is marked with a "+".

Below, two of the examples shown above are explicated:

Example 1 shows the wild type banding pattern. All wild type probes but none of the mutation probes display a signal; hence, the evaluation chart shows "+" in the three wild type columns and "-" in the three mutation columns. Accordingly, the boxes "FLQ-sensitive", "AP/CP-sensitive" and "EMB-sensitive" are marked with a "+".

In example 5, one of the *gyrA* wild type probes and the *embB* wild type probe are missing. Hence, the boxes "*gyrA* WT" and "*embB* WT" are marked with a "-". As none of the mutation probes are developed, these boxes are also marked with a "-". The *rrs* gene shows the wild type pattern. The strain is evaluated as FLQ- and EMB-resistant and AG/CP-sensitive.

Limitations

As with other diagnostic assays, the results of this test may only be interpreted in combination with additional laboratory and clinical data available to the responsible physician.

Use of this assay is limited to qualified personnel well trained in the test procedure and familiar with molecular biological methods.

The statements given in this manual concerning the capacities of the method refer to culture samples and pulmonary smear-positive clinical specimens.

The available data about extrapulmonary smear-positive clinical specimens is not sufficient to draw a conclusion about their applicability.

The presence of multiple bacterial species in the sample to be analyzed might hamper the interpretation of the test.

Prior to amplification, DNA has to be extracted from cultured bacteria or pulmonary smear-positive patient specimens using a suitable method. It must be ensured that the template DNA is efficiently amplified during the amplification reaction.

Performance evaluation of this assay was carried out using the HotStarTaq polymerase from Qiagen (Hilden, Germany). For DNA isolation from culture samples and clinical specimens the quick protocol described in chapter DNA Extraction was applied. Since the performance characteristics of this assay have not been validated for all polymerases and DNA isolation kits commercially available, the user is in charge of validating the applicability of polymerases or isolation methods other than those mentioned above.

As with any DNA-based assay, this test only screens the nucleic acid sequence and not the amino acid sequence. Therefore, it is possible that mutations that do not cause an amino acid exchange (silent mutations) will still produce the absence of one of the wild type probes.

Theoretically, a resistance can exist in spite of a wild type pattern. If, at investigation, the sample contains a strain that has developed a heteroresistance and the resistance of which is due to a mutation not covered by the mutation probes, the wild type pattern will appear. Similarly, if the sample contains

more than one *M. tuberculosis* complex strain (due to mixed culture or contamination) and one of these harbors a mutation not covered by the mutation probes, the wild type pattern will appear also.

The **GenoType® MTBDRsl** test only detects those resistances of the *M. tuberculosis* complex that have their origins in the *gyrA*, *rrs* and *embB* gene regions examined here. Resistances originating from mutations of other genes or gene regions as well as other resistance mechanisms against fluoroquinolones, aminoglycosides/cyclic peptides or ethambutol will not be detected by this test.

Resistance-mediating mutations within the *gyrA*, *rrs* and *embB* gene regions examined here which, due to their rareness, were not accessible for validation via patient samples, were only "detected" *in silico*. An *in silico* detection, however, does not exclude the possibility that some of these mutations cannot be detected *in vitro*.

The data given in table 2 regarding the cross-resistances between capreomycin, viomycin, kanamycin, and amikacin which vary depending on the mutations present, reflect the state of knowledge of Hain Lifescience and are based on the literature cited above. Recent data indicate that multiple mutations within and outside the sequences investigated with this assay may boost the intensity of phenotypic resistance or may shift the cross-resistances described above. Please note that effects due to multiple mutations outside the investigated sequences cannot be detected by this test.

The user must have or acquire information about the local mutation distribution pattern of the genes investigated with this test.

The test only works within the limits of the genomic regions the probes were chosen from. Potential sequence analysis remains to further investigations.

As with any detection system based on hybridization, the test system on hand bears the possibility that sequence variations in the genomic regions the primers and probes were chosen from but the detection of which the test was not designed for may lead to false results. Due to the high variability of bacterial genomes, it is possible that certain sub-types might not be detected. The test reflects the state of knowledge of Hain Lifescience.

Performance data of the assay can be requested through: www.hain-lifescience.com

Troubleshooting

Overall weak or no signals (including Conjugate Control zone)

- Room temperature too low or reagents not equilibrated to room temperature.
- No or too little amount of CON-C and/or SUB-C used.

Weak or no signals except for Conjugate Control zone

- Quality and/or quantity of extracted DNA do not allow an efficient amplification. Check amplicon on a 2% agarose gel. In case no amplicon is visible, repeat DNA extraction and amplification. If necessary, try a different DNA extraction method (see chapter DNA Extraction).
- Incubation temperature too high.

No homogeneous staining

- Strips were not completely immersed during incubation steps.
- Tray was not shaken properly.

High background color

- CON-C and/or SUB-C used too concentrated.
- Washing steps were not performed with the necessary care.
- Wash solutions too cold.

Unexpected result

- Wrong incubation temperature.
- Hybridization Buffer and/or Stringent Wash Solution were not properly prewarmed or mixed.
- Contamination of extracted DNA and/or amplification agents with extracted and/or amplified DNA. In case amplification agents are contaminated a negative control sample also shows the respective banding pattern.
- Contamination of neighboring wells by spillage during addition of Hybridization Buffer.
- Depending on the amount of amplified DNA used and on the specific reaction conditions, a strong and fast color development may occur. In such cases, discontinue the substrate incubation as soon as the signals are clearly visible in order to prevent the development of cross-hybridizing bands.
- No pure culture as starting material or more than one mutation in the tested strain.
- Silent mutation in probe region (see chapter Limitations).

Material Required but not Provided

- Absorbent paper
- Adjustable pipettes for 10, 20, 200, and 1000 µl
- Calibrated thermometer
- Disposable gloves
- Disposable sterile pipette tips with filter
- DNA extraction reagents for amplification use as well as necessary equipment
- Graduated cylinder
- PCR tubes, DNase and RNase free
- Shaking platform/**TwinCubator®**
- Shaking water bath/**TwinCubator®**
- Thermal cycler
(heating rate: 3°C/sec, cooling rate: 2°C/sec, precision: +/-0.2°C)
- Thermostable DNA polymerase with buffer (recommendation: hot start enzyme, extension rate: 2-4 kb/min at 72°C, half-life: 10 min at 97°C, 60 min at 94°C, amplification efficiency: >10⁵ fold)
- Timer
- Tweezers
- Water (molecular biology grade)
- Water bath

Kit Contents

Supplied

Membrane strips coated with specific probes (STRIPS)	12	96
Primer Nucleotide Mix (PNM) contains specific primers, nucleotides, <1% Dimethyl Sulfoxide, dye	0.5 ml	4 ml
Denaturation Solution (DEN) ready to use contains <2% NaOH, dye	0.3 ml	2.4 ml
Hybridization Buffer (HYB) ready to use contains 8-10% anionic tenside, dye	20 ml	120 ml
Stringent Wash Solution (STR) ready to use contains >25% of a quaternary ammonium compound, <1% anionic tenside, dye	20 ml	120 ml
Rinse Solution (RIN) ready to use contains buffer, <1% NaCl, <1% anionic tenside	50 ml	360 ml
Conjugate Concentrate (CON-C) concentrate contains streptavidin-conjugated alkaline phosphatase, dye	0.2 ml	1.2 ml
Conjugate Buffer (CON-D) contains buffer, 1% blocking reagent, <1% NaCl	20 ml	120 ml
Substrate Concentrate (SUB-C) concentrate contains Dimethyl Sulfoxide, substrate solution	0.2 ml	1.2 ml
Substrate Buffer (SUB-D) contains buffer, <1% MgCl ₂ , <1% NaCl	20 ml	120 ml
tray, evaluation sheet	1 of each	4 of each
manual, template	1 of each	1 of each